

INSTRUCTIONS on digitalMLPA Data Demo Set

NXtec™ D024-A1 KaryoProfiler

1. Purpose

This digitalMLPA data demo set serves to provide an overview of possible results that can be obtained with NXtec™ D024-A1 KaryoProfiler. It contains data files of eight samples with different genomic profiles and five reference samples. This data demo set also allows you to familiarise yourself with Coffalyser digitalMLPA™ and the way digitalMLPA results are presented by the software. The Sample Results PDF reports that are included in this demo set offer you the possibility to have a look at the sample results without analysing the data.

2. Description of samples included

Table 1 – Samples included in the D024-A1 data demo set

DNA Sample	Barcode used	Expected Results *
Coriell NA22976	BP04-49	1p: Heterozygous deletion affecting the probes for <i>AGRN</i> and <i>AJAP1</i> . 3q: Heterozygous duplication affecting the probes for <i>ATP11B</i> , <i>MAP3K13</i> , <i>ST6GAL1</i> , <i>LPP</i> , <i>TPRG1</i> , <i>IL1RAP</i> , <i>MB21D2</i> , <i>ACAP2</i> and <i>DLG1</i> . 9: Heterozygous duplication affecting the probes for <i>KANK1</i> , <i>VLDLR</i> , <i>JAK2</i> , <i>GLDC</i> , <i>PTPRD</i> , <i>LURAP1L</i> , <i>FREM1</i> , <i>BNC2</i> , <i>ADAMTSL1</i> , <i>MLLT3</i> , <i>IZUMO3</i> , <i>LINGO2</i> , <i>NFX1</i> , <i>GALT</i> , <i>PAX5</i> , <i>PIP5K1B</i> , <i>TRPM3</i> , <i>ANXA1</i> , <i>RORB</i> , <i>VPS13A</i> , <i>FRMD3</i> , <i>DAPK1</i> , <i>NXNL2</i> , <i>SYK</i> , <i>BICD2</i> , <i>FANCC</i> , <i>TDRD7</i> , <i>TGFBR1</i> , <i>ALDOB</i> , <i>SLC44A1</i> , <i>RAD23B</i> , <i>PTPN3</i> , <i>MUSK</i> , <i>KIAA1958</i> , <i>TNC</i> , <i>TRIM32</i> , <i>BRINP1</i> , <i>B3GALT9</i> , <i>RABGAP1</i> , <i>RALGPS1</i> , <i>NUP188</i> , <i>ABL1</i> , <i>OLFM1</i> and <i>NOTCH1</i> .
REF-Male-1	BP04-51	Reference sample: normal ratios for all probes.
Coriell NA03563	BP04-52	3: Heterozygous duplication affecting the probes for <i>TMEM45A</i> , <i>ALCAM</i> , <i>DUBR</i> , <i>DZIP3</i> , <i>NECTIN3</i> , <i>BOC</i> , <i>ZBTB20</i> , <i>TUSC7</i> , <i>UPK1B</i> , <i>POLQ</i> , <i>PARP14</i> , <i>UMPS</i> , <i>PLXNA1</i> , <i>RPN1</i> , <i>COL6A6</i> , <i>TMEM108</i> , <i>CEP63</i> , <i>NCK1</i> , <i>PIK3CB</i> , <i>CLSTN2</i> , <i>ATR</i> , <i>ZIC1</i> , <i>HLTF</i> , <i>MED12L</i> , <i>RAP2B</i> , <i>MME</i> , <i>LEKR1</i> , <i>GFM1</i> , <i>PPM1L</i> , <i>SLITRK3</i> , <i>SERPINI2</i> , <i>MECOM</i> , <i>TNIK</i> , <i>NLGN1</i> , <i>NAALADL2</i> , <i>TBL1XR1</i> , <i>PIK3CA</i> , <i>FXR1</i> , <i>ATP11B</i> , <i>MAP3K13</i> , <i>ST6GAL1</i> , <i>LPP</i> , <i>TPRG1</i> , <i>IL1RAP</i> , <i>MB21D2</i> , <i>ACAP2</i> and <i>DLG1</i> . 9p: Heterozygous deletion affecting the probe for <i>KANK1</i> .
Coriell NA00782	BP04-53	4p: Heterozygous duplication affecting the probes for <i>ADGRL3</i> , <i>EPHA5</i> , <i>CENPC</i> , <i>RUFY3</i> , <i>AFM</i> , <i>PARM1</i> , <i>SHROOM3</i> , <i>BMP2K</i> , <i>PRKG2</i> , <i>MRPS18C</i> , <i>CDS1</i> , <i>HERC3</i> , <i>CCSER1</i> , <i>GRID2</i> , <i>BMPR1B</i> , <i>ADH4</i> , <i>PPP3CA</i> , <i>NFKB1</i> , <i>DKK2</i> , <i>SEC24B</i> , <i>ENPEP</i> , <i>ANK2</i> , <i>UGT8</i> , <i>TRAM1L1</i> , <i>MYOZ2</i> , <i>PRDM5</i> , <i>SPRY1</i> , <i>FAT4</i> , <i>JADE1</i> , <i>PCDH10</i> , <i>RAB33B</i> , <i>TBC1D9</i> , <i>GAB1</i> and <i>HHIP</i> .
REF-Male-2	BP04-54	Reference sample: normal ratios for all probes.
Coriell NA10313	BP04-55	4q: Heterozygous duplication affecting the probes for <i>SH3D19</i> , <i>TRIM2</i> , <i>LRAT</i> , <i>PDGFC</i> , <i>FNIP2</i> , <i>NPY5R</i> , <i>KLHL2</i> , <i>SPOCK3</i> , <i>PALLD</i> , <i>HPGD</i> , <i>TENM3</i> , <i>CASP3</i> and <i>FAT1</i> . 7q: Heterozygous deletion affecting the probes for <i>RNF32</i> and <i>VIPR2</i> .
Coriell NA16445	BP04-57	17q: Heterozygous duplication affecting the probes for <i>MGAT5B</i> , <i>LGALS3BP</i> , <i>RPTOR</i> , <i>FN3KRP</i> and <i>TBCD</i> . 18p: Heterozygous deletion affecting the probe for <i>ENOSF1</i> .
Coriell NA07891	BP04-58	12q: Heterozygous duplication affecting the probes for <i>TMEM132D</i> , <i>ADGRD1</i> and <i>ZNF10</i> . 18q: Heterozygous deletion affecting the probes for <i>ELAC1</i> , <i>DCC</i> , <i>TCF4</i> , <i>WDR7</i> , <i>MALT1</i> , <i>CDH20</i> , <i>KDSR</i> , <i>CDH19</i> , <i>TMX3</i> , <i>GTSCR1</i> , <i>CBLN2</i> , <i>ZNF407</i> , <i>MBP</i> and <i>NFATC1</i> .
REF-Male-3	BP04-59	Reference sample: normal ratios for all probes.
Coriell NA09687	BP04-60	16p: Heterozygous deletion affecting the probe for <i>MSLN</i> .

		16q: Heterozygous duplication affecting the probes for <i>CALB2</i> , <i>ZFH3</i> , <i>RWD3</i> , <i>ADAMTS18</i> , <i>BCO1</i> , <i>CDH13</i> , <i>CRISPLD2</i> , <i>ZCCHC14</i> and <i>GALNS</i> .
<i>REF-Male-4</i>	BP04-62	<i>Reference sample: normal ratios for all probes.</i>
		4p: Heterozygous deletion affecting the probes for <i>PDE6B</i> , <i>ADD1</i> , <i>CYTL1</i> , <i>TBC1D14</i> , <i>GPR78</i> , <i>CLNK</i> , <i>CPEB2</i> , <i>LDB2</i> and <i>KCNIP4</i> .
Coriell NA03435	BP04-63	11p: Heterozygous duplication affecting the probes for <i>DEAF1</i> and <i>MMP26</i> .
<i>REF-Male-5</i>	BP04-64	<i>Reference sample: normal ratios for all probes.</i>

* Single probe aberrations may not be listed in the expected results.

3. Analysis in Coffalyser digitalMLPA software

The FASTQ file *FASTQ digitalMLPA Demo Set D024-A1.fastq* is generated by the Illumina® NGS platform. You can use this file to analyse the samples following the procedure in chapter **3.2 Analyse Your Data** as described in the [Coffalyser digitalMLPA User Manual](#).

Specific information has to be entered into Coffalyser digitalMLPA in several steps in section 3.2.1.1 *Configure general settings*. This concerns the steps listed below:

- Step 2: Leave the barcode lot on *From lot 03-009-yymmdd onwards (Default)*.
- Step 4a: Select the barcodes listed in Table 1.
- Step 7: Select *D024* from the product drop-down menu.
- Step 8: Select *A1-0125* from the sheet drop-down menu.
- Step 10: Select *Reference* from the drop-down menu in the column *type* for the reference samples (see Table 1). All other samples can be defined as *Test*.
- Step 11: Select *Pooled DNA source* from the drop-down menu in the column *options* for the reference samples. Leave this on *Default* for all other samples.